

Pediatric Patient Care Simulator

M56B

Product Supervision
 Medical Association of Osaka Prefecture
 Pediatric Medical Care Exploratory Committee



Pediatric Patient Care Simulator offers training of day-to-day pediatric patient care for healthcare professionals, teachers and others who need to acquire such skills.

Features

1. Life-like model of a 6 year old.
2. The simulator is able to sit in an upright position without support.
3. Oral, nasal and endotracheal airway suction
4. Both female and male genitalia are included for catheterization training

Training Skills

- Positioning
- Patient handling
- Airway suction:
 - oral, nasal or through a tracheostomy opening
- Percutaneous Endoscopic Gastrostomy (PEG) care
- Feeding tube management
- Male & Female catheterization



Airway Suction (oral)



Airway Suction (nasal)



Airway Suction (tracheostomy)



Feeding Tube Management



PEG care



Positioning, Patient Handling



Male & Female Catheterization



Anatomical Understanding

Set Includes

- 1 pediatric manikin
- 1 head & neck section model
- 1 male & female genitalia
- 1 syringe

- 1 water bag for catheter
- 1 strap-on IV venipuncture trainer
- 1 gym uniform
- 1 set of catheter (sample)

Manikin Size

125 cm

Replacement Parts

11355-030
 1 valve for catheterization