

SIMULATOR M60/72 Bathing and In-Patient Care Simulators "Fukutaro" and "Ofukusan"

M60.72 Bathing and In-Patient Care Simulators "Fukutaro" and "Ofukusan"



This model facilitates the various procedures of elderly care at anytime.

Because of the specially made skin of this model, bathing procedures are easy and harmless to the model.

Features:

1. Not only is the body seamless, but also waterproof to allow for repeated bathing of the model.
2. The hands and feet are movable in all directions, giving students accurate practice in the correct positioning and movement of patients.
3. This model is lightweight, providing for easy movement and storage.
4. The jaw of the male model is movable, so that oral care can easily be shown and practiced; furthermore, the hair of both models may be freely shampooed, combed, brushed, and blow-dried.



Practicable Techniques:

Denture and oral care practice (male model only), bathing and sponge bathing practice, rest, sleep, wheelchair, and bath positioning practice, massage, and others.

Specifications:

Male "Fukutaro" Dimensions: Approx. Length – 168cm (5'7")
Approx. Weight – 14kg (31lbs) 11274-000
Female "Ofukusan" Dimensions: Approx Length – 158cm (5'3")
Approx. Weight – 13kg (29lbs) 11287-000
Made with a metal skeleton frame and special resin skin.

 **KYOTO KAGAKU co.,LTD**

Call +81-75-605-2510 Fax +81-75-605-2519

Call +81-75-605-2520 Fax +81-75-605-2529

E-mail: kyoiku@kyotokagaku.co.jp

<http://www.kyotokagaku.co.jp>

If you need further information, please feel free to contact us at any time!

The specifications and looks may change without notice.