



- **Size**

12.5 x 13.6 x 7.8 (in)

- **Set Includes**

1 main body
1 stand for Hochstetter method
1 controller
1 power AC adapter
1 set of sensor pads

This simulator is designed to practice in the three ways of injection into the buttocks. It is made of a special resin that allows this model to be repeatedly injected without leaving a trace. The equipped indicator shows the accuracy of needle tip positioning and depth.

Features

1. The model is made from three individual layers : which are skin, muscular and skeleton. It allows identifying the injection spot.
2. The model is designed for three methods: dorsogluteal site, Hochstetter' s point and Clark 's point.
3. Attached indicator shows the accuracy of injection. If the injection spot is correct, the syringe can inject liquid.
4. It allows palpating the landmarks such as greater trochanter, superior anterior iliac spine, posterior superior iliac spine and sacrum.
5. Both skin and the muscle are replaceable.

Training Skills

- Gluteus maximus injection

- **Replacement Parts**

11213-010
skin



11213-020
sensor for quarter-trichotomy (pair)



11213-030
sensor for hochstetter method (pair)



11213-040
sensor for clark method (pair)



11213-050
controller

