

# Airway Suction Trainer

M85



**Features** A representation of human respiratory organs, designed for temporary airway suction and catheter insertion practice. Effective tool with life-like tactile qualities for training in airway suction for choking patients. Three catheter routes are prepared: oral, pharynx and tracheostomy opening. A sputum-like substance can be collected by proper suction. The airway and movement of the catheter tip are directly observed through the transparent side window.

1. Realistic anatomy includes the mouth, tongue, oral and nasal pharynx, epiglottis, trachea and esophagus.
2. Insertion of a catheter orally, nasally or through the tracheostomy opening.
3. Life-size organs facilitate learning proper catheter insertion length.
4. An observation window allows the passage of the catheter to be observed.
6. Realistic sticky simulated sputum.
7. The model is washable and can be dismantled easily for cleansing.

## Specifications

### Set includes:

- 1 head & neck model with airway anatomy
- 1 bottle of lubricant
- 2 catheters
- 1 bottle of simulated sputum
- 1 metal tray
- 1 syringe
- 1 storage case
- Model size: 45 x 23 x 15H cm
- Packing size: 61 x 31 x 43H cm, 7 kg
- Specifications are subject to change.

### Replacement parts

- 11229-020 replacement skin with trachea
- 11229-010 simulated sputum

