




* A wheelchair is not included.

Simple but high-quality, KEIKO reaches a wider range of trainees, including home-care nurses and caregivers

Keiko provides quality training for trainees and professionals who support long term care of patients in community and hospital.

Training Skills



- | | | | | |
|---|---|---|---|---|
| <p>Feeding</p> <ul style="list-style-type: none"> ● Tube Feeding (oral/nasal) <p>Elimination</p> <ul style="list-style-type: none"> ● Enema ● Urinary catheterization ● Stoma care ● Bedpan practice | <p>Exercise</p> <ul style="list-style-type: none"> ● Body positioning ● Wheelchair transfer <p>Hygiene</p> <ul style="list-style-type: none"> ● Oral care ● Bed bathing ● Partial bathing ● Change of clothes ● Perineal care ● Hair care* ● Grooming and appearance | <p>Respiratory Care</p> <ul style="list-style-type: none"> ● Oxygen inhalation ● Airway suction ● Postural drainage ● Tracheotomy care <p>Wound Care</p> <ul style="list-style-type: none"> ● Pressure sore care | <p>Drug Administration</p> <ul style="list-style-type: none"> ● Suppository | <p>Assessment</p> <ul style="list-style-type: none"> ● Procedure of vital signs taking <p>Pain Management</p> <ul style="list-style-type: none"> ● Fomentation <p>Others</p> <ul style="list-style-type: none"> ● Vaginal douche ● Postmortem care |
|---|---|---|---|---|

* optional

Features

BASIC PATIENT CARE

Oral Care

- Oral care on the bed and wheel chair
- Special care for the patient with feeding tube



- Brushing support
- Oral massage
- Denture care

Pressure sores care

- Understanding the levels of pressure sores
- Support to the patients having pressure sores



- (Stage 2) Scapula
- (Stage 3) Great trochanter
- (Stage 4) Sacrum bone
- (Unable to judge) Heel

Perineal care

- Anatomically correct genitals for detailed care



Body positioning and transfer

- Communication and appropriate handling to change the body positioning



- Changing body positioning
- Supporting to transfer the patient to a wheel chair

MEDICAL UNDERSTANDING

Airway Suction

- Understanding the procedure



- Airway suction: oral, nasal and tracheotomy (*tube insertion only)
- Tube feeding: oral, nasal and PEG (*tube insertion only)

Daily care with medical understandings



- Stoma care
- Bed bathing for the patient with central venous port
- Perineal care for the patient with urinal catheter



Specifications :

Set Includes:		Manikin Size:	Replacement Parts:	Optional parts:
1 manikin	1 male genitalia	149 cm, 15 kg	(11405-010) female genitalia	(12746-020) tracheotomy tube
1 denture	1 under pants		(11405-020) male genitalia	(11251-320)
1 wig	1 wrench		(11251-020) denture	simulated suppository set
1 oral cover for suctioning	1 instruction manual		(11229-060)	(11251-330) simulated feces set
1 female genitalia			oral cover for suctioning	(11405-030) pajama
				(11404-030) washable wig (elderly)