



### Perfect for a large number of trainees for IV

Basic IV Arm is realistic with fine detail, inexpensive and durable, perfect for large number of trainees to practice skills required in their courses.

Lightweight and portable, easy set up so that each trainee can have their own arm to practice.

The vein tube can easily be filled with simulated blood using a soft plastic injection bottle.

### Features

1. Flexible skin and palpable veins  
(replicates venous dilatation)
2. Realistic needle tip feeling and blood flashback
3. Injection pad is easily replaceable and durable.
4. No needle marks remain.
5. Hemorrhage can be stopped by compressing with the finger.
6. Adhesive plaster sticks well to the arm and can be easily removed.
7. Arm and supporter stand integrated so that users can set up easily.
8. Compact, lightweight and portable

### Training Skills

- Tourniquet application
- Venipuncture (including puncture with IV cannula)

### Specifications :

Set Includes:	Manikin Size:
4 arm model with supporter stand	W54 x D11 x D11 cm (each)
4 injection pad	
4 syringe (50ml)	
1 jar of coloring powder (red)	
1 spoon	
4 injection bottle for simulated blood	
1 carrying bag	