

Newborn Vital Signs Examination Simulator

Product No. **M58** 

An audible, palpable and thermometric simulator replays vital signs of a newborn child, now get CE marked. Heart sounds, arterial pulsation at the wrists, neck and temples, and also rectal temperature can be checked. Heart sounds are recorded from actual patient and are synchronized with arterial pulses. All vital signs can be controlled.

Production Supervision

Keiko Takahashi, Professor, College of Nursing Science,
Fukui Prefecture University



Auscultation




Palpation of pulses



Temperature taking

Training items

- Heart sounds
- Respiration
- Pulses:
 - superficial temporal,
 - carotid and radial arteries (bilateral)
- Body temperature (rectal thermometry)



Volume, heat and respiration rate and body temperature are controllable

Set includes:

- 1 newborn manikin
 - 1 control box
 - 1 thermometer
 - 1 baby clothes
 - 1 carrying case
- manikin size: Approx. 50 cm
- Specifications are subject to change.



www.kyotokagaku.com rw-kyoto@kyotokagaku.co.jp